



With CO2 connection



EcoFlushCap HybriX single use

Higher patient protection combined with lower costs



The innovative and single use **EcoFlushCap HybriX** supports the high demands on patient safety and significantly reduces the risk of cross-contamination. Connected to only one sterile water bottle system, the hybrid hose ensures safe flushing through the working channel and gentle flushing of the optics in the endoscope. The sterile, single use system is suitable for connection to all common flushing pump units through intelligent twin tubing. **EcoFlushCap HybriX** also supports CO2 insufflation thanks to the optional single use connection to CO2 insufflation devices.

Features

- ✔ Reduced risk of cross-contamination by eliminating the need to clean and reprocess reusable sterile water bottles
- ✔ The sterile system can be used for 24 hours
- ✔ For connection to all common flushing pump units
- ✔ With CO2 hose system
- ✔ With hose set and screw cap for all common sterile water bottles
- ✔ Flushing including check valve
- ✔ For Olympus endoscopes



Olympus/Fujifilm/Pentax



Lens cleaning

Ensures a clear view and optimal visibility

Olympus/Fujifilm/Pentax



External irrigation

Keeps the working channel free from residues

Olympus/Fujifilm/Pentax



CO2 insufflation

Reduces pressure and tension in patients

Twin tubing compatible with all common flushing pump units



Pieces per box	Manufacturer Product Code			
	For Olympus	For Fujifilm Series 500 / 600	For Fujifilm Series 700	For Pentax
	WITH 1.1 M CO2 CONNECTING HOSE			
15	WB-HYB-OL4	WB-HYB-FJ4	WB-HYB-FJ5	WB-HYB-PX6