

<b>Click`aV® Ligating Clips</b>				
Size/dimension	M	ML	L	XL
Stability type on the tissue	Crush clips compression exerted by the clip legs on the tissue	Crush clips compression exerted by the clip legs on the tissue	Crush clips compression exerted by the clip legs on the tissue	Crush clips compression exerted by the clip legs on the tissue
Distance between clip arms [mm]	8.7	11.9	16	19.4
Clip thickness [mm]	0.9	0.85	1.3	1.3
Recommended ligated structure size [mm]	2 - 7	3 - 10	5 - 13	7 - 16
<b>Click`aV Plus™ Ligating Clips</b>				
Size/dimension	ML	L	XL	XXL
Stability type on the tissue	Anchoring of the clip in the tissue	Anchoring of the clip in the tissue	Anchoring of the clip in the tissue	Anchoring of the clip in the tissue
Distance between clip arms [mm]	12	16	19	28.1
Clip thickness [mm]	0.85	1.3	1.3	2
Recommended ligated structure size [mm]	3 - 10	5 - 13	7 - 16	10 - 22