


Irrigation sets for urology/gynaecology

U91 (Spike B86, MLL)



1. The special design of the spike enables easy penetration of container ports and offers a safe connection to the container. Non-vented version. The protection cap maintains sterility of the set.
2. Roller clamp with on/off function.
3. Transparent flexible PVC tubing with plasticiser TOTM 4,8/7,2 mm (ID/OD) (large bore).
4. Male Luer-Lock adapter guarantees a safe connection to the device. The protection cap maintains the sterility of the set.



Disclaimer:

The above is indicated only for information. This document is not a specification. The values given are subject to the usual fluctuations. Please follow the instructions for use by the medication manufacturer and the actual SPC (Summary of Product Characteristics). Common hygiene measures and working methods are to be followed at all times.

Subject to alterations without prior notice!

REF	Product description	Packaging	GTIN	Filling/Storage volume
76.5705	Spike B86, 200 cm PVC (TOTM) tubing 4,8/7,2 mm (ID/OD), MLL adapter, FlowStop	1 30 960	04055698023929 04055698023936 -	